








**Laryngoscopes**  
**Laringoscópios**  
**Laryngoscopes**

Tel# 0092-52-4260785  
Fax# 0092-52-4266804  
Email. [Info@benison.com.pk](mailto:Info@benison.com.pk)

**Tongue Holding Forceps, Catheter Introducing Forceps**



**Mouth Gags**



### Laryngoscope Sets, McIntosh



31.23.40 – 31.23.55



31.43.40 – 31.43.55



31.22.00 – 31.23.05  
31.42.00 – 31.45.05

#### MCINTOSH (CONVENTIONAL)

| Wooden Box | Artificial Leather<br>Emergency Bag | Contents                                |
|------------|-------------------------------------|---|
| 31.22.00   | 31.23.40                            | Set of 3 blades No. 2,3,4               |
| 31.22.03   | 31.23.53                            | Set of 4 blades No. 0,1,2,3 (Pediatric) |
| 31.23.04   | 31.23.54                            | Set of 4 blades No. 1,2,3,4             |
| 31.23.05   | 31.23.55                            | Set of 5 blades, No. 1,2,3,4,5          |

All sets include 1 standard battery handle  
and 2 spare lamps

#### MCINTOSH (FIBEROPTIC)

| Wooden Box | Artificial Leather<br>Emergency Bag | Contents                                |
|------------|-------------------------------------|---|
| 31.42.00   | 31.43.40                            | Set of 3 blades No. 2,3,4               |
| 31.42.03   | 31.43.50                            | Set of 4 blades No. 0,1,2,3 (Pediatric) |
| 31.43.04   | 31.43.54                            | Set of 4 blades No. 1,2,3,4             |
| 31.43.05   | 31.43.55                            | Set of 5 blades, No. 1,2,3,4,5          |

All sets include 1 standard battery handle  
and 2 spare lamps

**Laryngoscope Sets, Wis-Forregger, Miller**



31.26.20 – 31.26.26  
31.46.40 – 31.46.56

**WIS-FOREGGER (CONVENTIONAL)**

| Wooden Box | Artificial Leather<br>Emergency Bag | Contents                                |
|------------|-------------------------------------|---|
| 31.26.00   | 31.26.20                            | Set of 3 blades No. 2,3,4               |
| 31.26.04   | 31.26.24                            | Set of 4 blades No. 0,1,2,3 (Pediatric) |
| 31.26.05   | 31.26.25                            | Set of 4 blades No. 1,2,3,4             |
| 31.26.06   | 31.26.26                            | Set of 5 blades, No. 1,2,3,4,5          |

All sets include 1 standard battery handle  
and 2 spare lamps

**WIS-FOREGGER (FIBRE OPTIC)**

| Wooden Box | Artificial Leather<br>Emergency Bag | Contents                                |
|------------|-------------------------------------|---|
| 31.46.00   | 31.46.40                            | Set of 3 blades No. 2,3,4               |
| 31.46.04   | 31.46.54                            | Set of 4 blades No. 0,1,2,3 (Pediatric) |
| 31.46.05   | 31.46.55                            | Set of 4 blades No. 1,2,3,4             |
| 31.46.06   | 31.46.56                            | Set of 5 blades, No. 1,2,3,4,5          |

All sets include 1 standard battery handle  
and 2 spare lamps



31.28.20 – 31.28.28  
31.48.40 – 31.48.56

**MILLER (CONVENTIONAL)**

| Wooden Box | Artificial Leather<br>Emergency Bag | Contents                                |
|------------|-------------------------------------|---|
| 31.28.00   | 31.28.20                            | Set of 3 blades No. 1,2,3               |
| 31.28.04   | 31.28.24                            | Set of 4 blades No. 0,1,2,3 (Pediatric) |
| 31.28.05   | 31.28.25                            | Set of 4 blades No. 1,2,3,4             |
| 31.28.06   | 31.28.28                            | Set of 5 blades, No. 1,2,3,4,5          |


All sets include 1 standard battery handle  
and 2 spare lamps

**MILLER (FIBRE OPTIC)**

| Wooden Box | Artificial Leather<br>Emergency Bag | Contents                                |
|------------|-------------------------------------|---|
| 31.48.00   | 31.48.40                            | Set of 3 blades No. 1,2,3               |
| 31.48.04   | 31.48.54                            | Set of 4 blades No. 0,1,2,3 (Pediatric) |
| 31.48.05   | 31.48.55                            | Set of 4 blades No. 1,2,3,4             |
| 31.48.06   | 31.48.56                            | Set of 5 blades, No. 1,2,3,4,5          |

All sets include 1 standard battery handle  
and 2 spare lamps

**Mcintosh Blades, Conventional**



**31.24.10**  
Fig. 0, 55 mm  
For Premature Babies



**31.24.11**  
Fig. 1, 70 mm  
For Infants



**31.24.12**  
Fig. 2, 90 mm  
For Children



**31.24.13**  
Fig. 3, 110 mm  
For Adults, Medium




**31.24.14**  
Fig. 4, 130 mm  
For Adults, Large



**31.24.15**  
Fig. 5, 150 mm  
For Adults, Extra Large

Mcintosh Blades, Fibre Optic cold light



**31.44.10**  
Fig. 0, 55 mm  
For Premature Babies



**31.44.11**  
Fig. 1, 70 mm  
For Infants



**31.44.12**  
Fig. 2, 90 mm  
For Children



**31.44.13**  
Fig. 3, 110 mm  
For Adults, Medium



**31.44.14**  
Fig. 4, 130 mm  
For Adults, Large

Wis-Foregger Blades, Conventional



**31.26.10**  
Fig. 0, 55 mm  
For Premature Babies



**31.26.11**  
Fig. 1, 70 mm  
For Infants



**31.26.12**  
Fig. 2, 90 mm  
For Children



**31.26.13**  
Fig. 3, 110 mm  
For Adults, Medium



**31.26.14**  
Fig. 4, 135 mm  
For Adults, Large



**Wis-Foregger Blades, Fibre Optic - Cold light**



**31.46.10**  
Fig. 0, 55 mm  
For Premature Babies



**31.46.11**  
Fig. 1, 70 mm  
For Infants



**31.46.12**  
Fig. 2, 90 mm  
For Children



**31.46.13**  
Fig. 3, 110 mm  
For Adults, Medium



**31.46.14**  
Fig. 4, 135 mm  
For Adults, Large

**Miller Blades, Conventional**



**31.28.30**  
Fig. 0, 55 mm  
For Premature Babies



**31.28.31**  
Fig. 1, 80 mm  
For Infants



**31.28.32**  
Fig. 2, 130 mm  
For Children



**31.28.33**  
Fig. 3, 170 mm  
For Adults, Medium



**31.28.34**  
Fig. 4, 180 mm  
For Adults, Large

Miller Blades, Fibre Optic- Cold Light



**31.48.30**  
Fig. 0, 55 mm  
For Premature Babies



**31.48.31**  
Fig. 1, 80 mm  
For Infants



**31.48.32**  
Fig. 2, 130 mm  
For Children



**31.48.33**  
Fig. 3, 170 mm  
For Adults, Medium



**31.48.34**  
Fig. 4, 180 mm  
For Adults, Large

**Battery Handles & Lamps**



FIBRE OPTIC  
**31.30.40**  
ø 24 mm  
11 cm / 4 1/4"



STANDARD  
**31.31.11**  
ø 19 mm  
13 cm / 5"



STANDARD  
**31.32.12**  
ø 28mm  
14 cm / 5 1/2"



FIBRE OPTIC  
**31.32.62**  
ø 28mm  
14 cm / 5 1/2"